


IS THE HAIRMAX RIGHT FOR YOU?

HairMax® is most effective in women and men in early to moderate stages of hair loss and thinning. If your hair loss looks like the images below, HairMax is recommended for you.

FEMALE



MALE



DOCTOR RECOMMENDED

Physicians worldwide recommend HairMax Laser Light Therapy to their patients. HairMax is supported by a Medical Advisory Board of over 30 established key authoritative figures in the fields of dermatology and hair restoration.

FREQUENTLY ASKED QUESTIONS

How often do I need to use my LaserBand 82?

It is recommended to use your LaserBand three times a week for 90 seconds, with at least a day between each treatment.

Will results improve if I use my LaserBand 82 more than 3 times a week?

Three times a week has been found to be the optimal treatment time. Using your LaserBand more than 3 times a week has not shown any additional benefits.

How long before I see results?

Most users start to see benefits in as little as 16 weeks.

How long will I need to continue hair loss treatments with the LaserBand 82?

To keep seeing benefits use your LaserBand as directed on an ongoing basis. If you do not use the device as recommended, there is a strong possibility that you will not see benefits or your hair loss will regress to where you were prior to using the LaserBand.

Can I use my LaserBand 82 on wet or dry hair?

You can use your device on dry or damp hair. It is important that your scalp be clean and free of styling products. A build-up of topicals such as oils, gels, sprays, etc. can block the laser light from effectively reaching the hair follicles.

Do I have to use any special shampoos or products with my LaserBand 82?

We recommend using the HairMax for den•si•ty collection. These sulfate-free, low pH products are gentle on your thinning hair and complement the laser treatment.




Lexington International, LLC  
777 Yamato Road, Suite 105, Boca Raton, FL 33431  
Tel. +1.561.417.0200 | U.S. Toll Free 1.800.9REGROW  
FAX. +1.561.892.0747

[www.HairMax.com](http://www.HairMax.com)

LIT-CONBROLEX-LB82-00 08/15 © Copyright 2015 Lexington Intl., LLC. All rights reserved.



THE NEXT  
BREAKTHROUGH IN  
HAIR GROWTH




HairMax®  
LaserBand 82

DESIGN BY *pininfarina*

**HAIRMAX®**  
The science of hair growth.™

HOW LASER LIGHT THERAPY  
PROMOTES HAIR GROWTH

HairMax® works through a process called Photobiostimulation, similar to photosynthesis in plants.




NOURISH. STIMULATE. GROW.

LIGHT IS ENERGY.

HairMax delivers therapeutic, nourishing light energy to your hair follicles to activate them for stronger growth. This boost of light revitalizes and enlivens weakened hair follicles to stimulate growth factors and reverse the thinning process. This will help extend and restore your hair's natural growth cycle, resulting in hair growth with increased density, fullness and vibrancy.



# New Hair Growth Beyond What You Thought Possible!



82 Medical Grade Lasers

DESIGN BY *pininfarina*

- Stimulates Hair Growth
- Treats Hair Loss
  - Reverses the Thinning Process
  - Increases Density and Fullness
  - Revitalizes Damaged Hair
  - No Harmful Side Effects
  - Drug-Free Treatment

**CLINICALLY  
TESTED**

**FDA  
CLEARED\***

**DOCTOR  
RECOMMENDED**

## CLINICALLY PROVEN

HairMax has been proven safe and effective in 7 clinical studies with a total of 460 participants and is the only laser hair loss device with clinical studies published in 5 medical peer reviewed journals. Participants experienced significant hair growth with an average increase in hair count of 129 additional hairs per square inch.

**90%**  
of participants in  
clinical studies  
experienced significant  
hair growth

\*HairMax laser devices are indicated to treat Androgenetic Alopecia, and promote hair growth in males who have Norwood Hamilton Classifications of IIa to V and in females who have Ludwig (Savin) I-4, II-1, II-2, or frontal patterns of hair loss and who both have Fitzpatrick Skin Types I to IV. Laser Light - Avoid direct eye exposure.

## HairMax® LaserBand 82

Stimulating hair growth has never been easier, faster or more sophisticated. Featuring 82 medical grade lasers, the HairMax LaserBand provides a **super-fast 90 second treatment**. The LaserBand's unique cradle design and exclusive hair parting teeth provide optimal delivery of nourishing laser light to energize your hair follicles and regrow your hair!

## SEEING IS BELIEVING

BEFORE



AFTER



**Female:** Significant new hair growth in front temple area.

BEFORE




AFTER



**Male:** Significant new hair growth, density and fullness.

## Design by Pininfarina

HairMax has collaborated with the leading Italian design firm Pininfarina, to create a refined and distinctive iconic design for the newest range of HairMax Laser devices. Combining form with function, technology with innovation, the HairMax-Pininfarina line of laser hair growth devices are exquisitely designed and unequalled in style, usability and most importantly, effectiveness.



NOURISH.  
STIMULATE.  
GROW.

